




KN95 Protective Face Mask



N95 Particulate Respirator Mask



3 Ply Single Use Medical Mask

Face Masks

Our range of CoShield face masks offer high grade protection to reduce your risk of contamination.

We offer a range of different mask sizes and specifications, from our simple single-use medical mask to our particulate respirator mask.

3 Ply Single Use Medical Mask

Features:

- Has synthetic blood penetration resistance
- 50 pieces per box
- Designed for single use
- Colour: blue/white
- Certificates & Standards
 - CE
 - Spec Sheet

KN95 Protective Face Mask

Features:

- Medical grade
- Ear loop style
- Colour: white
- Size 17.5 x 9.5cm
- Certificates:
 - CE
 - EN14683:2005
 - FDA
 - GTS
 - FFP2
 - EN149

N95 Particulate Respirator Mask

Features:

- Non-woven fabric
- Melt blown filter
- Size: 15 x 11cm
- Certificates:
 - FFP2
 - NIOSH
 - CE

Product Codes:

3ply: 03-04-14281

KN95: 03-04-14280

N95: 03-04-14230