



# OXYGEN CONCENTRATORS

2021 Global Catalog of Medical Oxygen Concentrator Equipment

# CoShield OXC Oxygen Concentrators




Medical oxygen supplied by the concentrator is beneficial to cure the respiratory disease or heart and blood vessel system, chronic pulmonary system, and other oxygen lacking symptoms.


## Key Benefits

1. Strengthens your heart, reducing the risk of heart attacks and help with cardiovascular activity.
2. Speeds up the body's recovery after physical exertion.
3. Stimulate brain activity, increase memory capacity.
4. Prevent disease, relieve hangovers.
5. Reduce stress and anxiety.
6. Alleviates tension headaches.
7. Remedy for irregular sleeping patterns.

## Interior



## Exterior



# OXC-005

## Product Technical Data



Low noise  
Takes up less space



USA imported humidifier  
6PSI pressure release



Innovative big  
LCD display

## Specifications

Flow Rate	0-5LPM
Purity	93% ( $\pm 3\%$ )
Outlet Pressure	0.04-0.07MPA(6-10PSI)
Sound Level	$\leq 38\text{db}$
Power Consumption	$\leq 350\text{w}$
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm High&low pressure alarm
Net Weight	15Kgs
Size	348×280×510 mm

## Options

1. Nebulizer (Atomization)
2. Low Purity alarm
3. Low purity alarm with low flow alarm
4. High temperature alarm
5. Pulse oximeter
6. Maintenance alarm
7. Reset Switch
8. Remote control monitoring system (WIFI/Bluetooth)

## Certifications

- ✓ FDA approved
- ✓ CE approved
- ✓ EN ISO13485:2016

# OXC-010


## Product Technical Data



**DUAL FLOWS**  
USA imported humidifier  
6PSI pressure release



Innovative "cooling"  
technology to improve  
reliability and life of the unit



**Big LCD display**  
24H\*365D work continuously

## Specifications

Flow Rate	0-10LPM
Purity	93% (±3%)
Outlet Pressure	0.04-0.07MPa(6-10PSI)
Sound Level	≤50db
Power Consumption	≤880w
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm High&low pressure alarm
Net Weight	27Kgs
Size	365×375×600 mm

## Options

1. Nebulizer (Atomization)
2. Low Purity alarm
3. Low purity alarm with low flow alarm
4. High temperature alarm
5. Pulse oximeter
6. Maintenance alarm
7. Reset Switch
8. Remote control monitoring system (WIFI/Bluetooth)

## Certifications

- ✓ FDA approved
- ✓ CE approved
- ✓ EN ISO13485:2016





[www.coshield.com](http://www.coshield.com)  
[hello@coshield.com](mailto:hello@coshield.com)

UNITED STATES CANADA CARIBBEAN ARGENTINA UNITED KINGDOM FRANCE THAILAND MALAYSIA CHINA AUSTRALIA NEW ZEALAND